

9301 Fast Preamplifier

- Used with photomultipliers, electron multipliers, etc., in photon or ion counting applications
- 1.5-ns rise time
- · Voltage gain of 10
- Output of ±0.7 V into 50 Ω
- · Compact and lightweight



The ORTEC Model 9301 low-noise, fast-rise-time Preamplifier has been designed for use with photomultipliers, electron multipliers, and other detectors employed in photon counting, ion counting, or fast-timing applications. When connected to the detector, the 50- Ω input resistance of the Model 9301 provides a load resistance for the detector output current pulse.

In addition to its fast rise time of 1.5 ns, this preamplifier has a voltage gain of 10 and an output impedance of $50~\Omega.$ Because of its compact size and light weight, the Model 9301 is ideal for mounting directly or close to a detector. Consequently, low-level signals, which would otherwise be susceptible to pickup of noise or interference, are boosted to a suitable level for cable connection to the main amplifier. Connection to a power supply is through the 3-m (10-ft) long captive power cable furnished with the Model 9301.

Model 9301 is fitted with a power cable connector (Amphenol 17-20090) that is compatible with other ORTEC NIMstandard modules.

Specifications

PERFORMANCE

INPUT IMPEDANCE 50 Ω .

VOLTAGE GAIN 10 (±2%) noninverting.

RISETIME <1.5 ns.

INPUT RMS NOISE EQUIVALENT

OUTPUT IMPEDANCE Typically 50 Ω . OUTPUT DYNAMIC RANGE >±0.7 V into 50 Ω .

NONLINEARITY <±1%.

TEMPERATURE GAIN INSTABILITY <±0.1%/°C.

INPUT CONNECTOR Front-panel BNC.

OUTPUT CONNECTOR Rear-panel BNC.

ELECTRICAL AND MECHANICAL

POWER REQUIRED +12 V, 30 mA; -12 V, 30 mA.

WEIGHT

Net 0.17 kg (6 oz).

Shipping 0.9 kg (~2 lb).

DIMENSIONS 3.1 X 5.0 X 7.3 cm (1.25 X 2 X 2.875 in.) plus 3-m (10-ft) cable.

Ordering Information

To order, specify:

Model Description
9301 Preamplifier

Suggested cable accessories:

C-25-1/2 RG-58A/U 50-Ω Cable with

two BNC male plugs; 6-in.

length

C-25-12 RG-58A/U 50-Ω Cable with

two BNC male plugs; 12-ft

length



www.ortec-online.com

Tel. (865) 482-4411 • Fax (865) 483-0396 • ortec.info@ametek.com 801 South Illinois Ave., Oak Ridge, TN 37831-0895 U.S.A. For International Office Locations, Visit Our Website



Specifications subject to change